



## User Manual

**Serum amyloid A proteins (SAA) (Rhesus macaque)**

**Cat. No. EEOPP-19012**

**Description:**

Serum amyloid A proteins (SAA) represents a family of apolipoproteins that circulates in association with high-density lipoproteins (HDL). The level of apo-SAA, normally 1-5  $\mu\text{g}/\text{mL}$  in plasma, increases 500-1000 fold within 24 hours of an inflammatory stimulus and, under these conditions, is the most abundant HDL apo-lipoprotein.

**Source:**

*Escherichia coli*

**Unit:**

10  $\mu\text{g}$

**Reconstitution:**

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Reconstitute in sterile distilled water or aqueous buffer containing 0.1% BSA to a concentration of 0.1-1.0  $\text{mg}/\text{mL}$ . Stock solutions should be apportioned into working aliquots and stored at  $< -20^{\circ}\text{C}$ . Further dilutions should be made in appropriate buffered solutions.

**Formulation:**

Lyophilized from a 0.2  $\mu\text{m}$  filtered concentrated solution in PBS, pH 7.4.

**Storage:**

This lyophilized preparation is stable at  $2-8^{\circ}\text{C}$ , but should be kept at  $-20^{\circ}\text{C}$  for long term storage, preferably desiccated. Upon reconstitution, the preparation is stable for up to one week at  $2-8^{\circ}\text{C}$ . For maximal stability, apportion the reconstituted preparation into working aliquots and store at  $-20^{\circ}\text{C}$  to  $-70^{\circ}\text{C}$ . Avoid repeated freeze/thaw cycles.

**Molecular Weight:**

Approximately 11.8 kDa, a single non-glycosylated polypeptide chain containing 104 amino acids.

**Endotoxin:**

Less than 1 EU/ $\mu$ g of SAA1 as determined by LAL method.

**Usage:**

This material is offered by Cyagen Biosciences for research, laboratory or further evaluation purposes. FOR RESEARCH USE ONLY. NOT INTENDED FOR ANY ANIMAL OR HUMAN THERAP EUTIC OR DIAGNOSTIC USE.

**Biological Activity:**

The biologically active determined by a chemotaxis bioassay using human monocytes is in a concentration range of 10-100 ng/mL.

**Physical Appearance:**

Sterile filtered white lyophilized (freeze-dried) powder.

**AA Sequence:**

RSWFSFLGEA YDGARDMWRA YSDMKEANYK NSDKYFHARG NYDAAQRGPG  
GVWAAEVISD ARENIQKLLG RGAEDTLADQ AANEWGRSGK DPNHFRPAGL PEKY

**Purity:**

> 97% by SDS-PAGE and HPLC analyses.

Material Safety Data Sheets (MSDSs) are available upon request.

The Certificate of Analysis (COA), which provides detailed quality control information for each product, is also available at the Cyagen website.

Cyagen Biosciences reserves all rights on the technical documents of its culture products. No part of this document may be reproduced or adapted for other purposes without written permission from Cyagen Biosciences.